




SPECTRONICS INDIA

AN ISO 9001:2015
CERTIFIED COMPANY

SI-165 DIGITAL COLONY COUNTER

Our diverse range of Digital Colony Counter offered by us is of premium quality and is fabricated using quality raw material. With the assistance of our designers and experts we are capable to offer customized solutions to our clients as per their specifications in stipulated time frame.



SALIENT FEATURES

-  Microprocessor Based counting control Memory against power failure. (Counted Data Not Lost)
-  Audible Confirmation of each count
-  110mm Dia Magnifying Glass
-  Uniform glare free illumination
-  Digital Readout (4 Digits) from 0 - 9999
-  Auto Marker pen
-  Wolffhuegel glass grid with focusing facility

ACCESSORIES

- Marking Pen
- Magnifier Lens
- Dust Cover
- Manual

TECHNICAL SPECIFICATIONS

PARAMETER	MODEL SI-165
DIGITAL DISPLAY	4 digit, 0 - 9999 Maximum Count
MEMORY	Microprocessor based with facility of Counting control memory against Power Failure
DISH SIZE	110 mm
MAGNIFICATION	X 1.7
DIMENSIONS	L 265 X B 230 X H 130 mm (Approx.)
WEIGHT	2 Kg. (Approx.)
POWER	230V ± 10% AC, 50Hz, 40W
ACCESSORIES	Marking Pen - (1 No.), Magnifier Lens - (1 No.)



info@spectronicsindia.in



Spectronicsindia.com



+91 9417350423

