



**SPECTRONICS INDIA**

AN ISO 9001:2015 &  
CE CERTIFIED COMPANY

# SI-194

## MICROPROCESSOR BASED DIGITAL CONDUCTIVITY METER

There is a wide range of Microprocessor based Digital Conductivity/Temp. Meter manufactured by us. These are made under the introspection of the experts which is why these are superior in quality as compared to the rest of the competitors. The competitive prices of these products make them more desirable in the market. The raw material used is of unmatched quality and provides for the better and prolonged usage of these products. The compact designs make them easy to operate and saves lot of time. Digital Conductivity Meters available with us assure you of durability accompanied with least maintenance.



### SALIENT FEATURES

-  HIGH RELIABILITY
-  BETTER PERFORMANCE
-  LED DISPLAY OF 3 1/2 DIGIT
-  FACILITY OF CELL CONSTANT ADJUSTMENT
-  AUTO RANGE
-  FACILITY OF AUTOMATIC TEMPERATURE COMPENSATION
-  FABRICATED USING QUALITY RAW MATERIAL
-  AVAILABLE AT OPTIMUM PRICES

### ACCESSORIES

- Conductivity Cell
- Temp. Sensor
- Operation Manual
- Cell Stand
- Dust Cover

SPECIFICATIONS	SI-194 TABLE MODEL
Model	SI-194 Table Model
Display	3 1/2 Digital LED
Auto Range	0-200µS 0-2,0-20,0-200 and 0-1000mS
Temperature Compensation	Manual 0-50°C and Automatic
Cell Constant Adjustment	Adjustable
Conductivity Cell	Platinum DIP Type
Resolution	0.1µS
Accuracy	± 1% FS ± 1 Digit
Power	230v ± 10% ac, 50 Hz
Dimensions	220 x 190 x 100 mm
Weight	1 kg. (Approx.)