



SPECTRONICS INDIA

ISO 9001:2015 & CE CERTIFIED COMPANY

DIGITAL MELTING POINT APPARATUS

Digital Melting Point Apparatus is a compact instrument designed for accurate determinations of Melting Points of solid samples.

MODEL SI-252



SALIENT FEATURES

- Cylindrical Silicon Oil Bath
- Built-in Magnetic Stirrer with electronic speed controller.
- Electronic Controller for adjusting the heating rate.
- Digital Display of Temperature with PT-100 Sensor.
- Glare-free back-ground light with adjustable light Intensity.
- Provision for holding the melting point display



Digital Readout



For Temperature
1°C



Silicon Oil-Bath



In-built
Magnetic Stirrer

TECHNICAL SPECIFICATION

Description	SI-252
Temperature Range	20 above RT-275°C
Readability	1°C
Display	4 Digit , 7 Segment LED Display
Heating Rate Control	Electronic
Stirrer Speed Control	Electronic
Temperature Sensor	PT-100
Weight	12 kg. Approx.

PRODUCT RANGE

Double Beam UV/VIS Spectrophotometer, Single Beam UV/VIS Spectrophotometer, Flame Photometer, Colorimeter, pH Meter, Conductivity Meter, DO Meter, Salinity Meter, Turbidity Meter, Nephelometer, Colony Counter, Auto Karl- Fischer, Friability Apparatus, Melting Point Apparatus, Temperature Indicator, Water & Soil Analysis Kit, Portable pH Meter, Student Microscope, Binocular Microscope, Medical Microscope and Portable Instruments.



gmail - info@spectronicsindia.in



website spectronicsindia.com