



# SPECTRONICS INDIA

ISO 9001:2015 & CE CERTIFIED COMPANY

## DIGITAL MELTING POINT APPARATUS

SI-253



Digital Readout



For Temperature 0.1 °C



Silicon Oil-Bath



In-built Magnetic Stirrer

### TECHNICAL SPECIFICATION

Description	SI-253
Temperature Range	20 above RT-2750C
Readability	0.10
Display	4 Digit , 7 Segment LED Display
Heating Rate Control	Electronic
Stirrer Speed Control	Electronic
Temperature Sensor	PT-100
Weight	12 kg. Approx.

### PRODUCT RANGE

Double Beam UV/VIS Spectrophotometer, Single Beam UV/VIS Spectrophotometer, Flame Photometer, Colorimeter, pH Meter, Conductivity Meter, DO Meter, Salinity Meter, Turbidity Meter, Nephelometer, Colony Counter, Auto Karl- Fischer, Friability Apparatus, Melting Point Apparatus, Temperature Indicator, Water & Soil Analysis Kit, Portable pH Meter, Student Microscope, Binocular Microscope, Medical Microscope and Portable Instruments.



gmail - [info@spectronicsindia.in](mailto:info@spectronicsindia.in)



website [spectronicsindia.com](http://spectronicsindia.com)