



SPECTRONICS INDIA

AN ISO 9001:2015 &
CE CERTIFIED COMPANY

SI-222

DIGITAL NEPHELO METER

We are offering our superior range of turbidity meter for our clients based globally. These are manufactured with the help of our designers and engineers using high technology and state-of-the-art facilities. Also, we are able to custom design the range in accordance with the specifications of the clients.



SALIENT FEATURES

- Highly Stable & Accurate
- 3 ½ Digit LED Display
- 30 mm flat bottom Test Tube
- 0 to 200 NTU

ACCESSORIES

- Test tubes
- Operation Manual
- Dust Cover

MODEL	SI-222
Type	Table Model
Display	3½ Digital LED
Detector	Photocell / Photodiode
Range	0 to 200 NTU
Resolution	0.1 NTU
Accuracy	±2% FS ± 1 Digit
Sample System	30 mm Clear Glass Test Tube
Light Source	6V Tungsten Lamp or Bright Blue LED
Power	220V ±10%, 50Hz AC
Dimensions	230mm x 180mm x 110mm (L x B x H)