



SPECTRONICS INDIA

AN ISO 9001:2015 &
CE CERTIFIED COMPANY

SI-139 DIGITAL pH METER

ACCURATE • RELIABLE • EASY TO USE

SALIENT FEATURES



**HIGHLY
STABLE &
ACCURATE**



**3 ½ DIGIT
LED DISPLAY**



**SLOPE
CONTROL**



TECHNICAL SPECIFICATIONS

Parameter	Range / Details	Specification
Model	139	Table Model
Display	3½ Digit LED	
Range	pH mV	0 to 14.00 0 to + 1999
Resolution	pH mV	0.01 1
Accuracy	pH mV	±0.01 ±1
Temperature	-	
Compensation	Manual 0 to 100°	
Slope Correction	80 to 120%	
Power	230V ±10%	
Dimensions	275 x 175 X 76	

ACCESSORIES



Combination pH
Electrode



Buffer Tablets
4pH & 7pH



Dust Cover



Stand for
Electrode



The diverse range of Digital pH meter catered by us is of top-notch quality and is manufactured using quality raw material and advanced technology. Our experts and engineers are capable enough to maintain pace with the latest trend of the marketplace. Further, we can custom design the range as per the specifications of the clients in the stipulated time frame.



spectronicsindia@gmail.com



spectronicsindia.in



+91 9417350423