



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

SI-143 TABLE MODEL

DELUXE pH METER



The diverse range of Deluxe pH meter catered by us is of top-notch quality and is manufactured using quality raw material and advanced technology. Our experts and engineers are capable enough to maintain pace with the latest trend of the marketplace. Further, we can custom design the range as per the specifications of the clients in the stipulated time frame.

SALIENT FEATURES

- Highly Stable & Accurate
- 3 ½ Digit LED Display
- ATC Facility available
- Slope Control

ACCESSORIES

- Combination pH Electrode,
- Buffer Tablets 4pH & 7pH,
- Dust cover and Stand for Electrode.



TECHNICAL SPECIFICATIONS

PARAMETER	RANGE	SI-143 TABLE MODEL
Display	–	3 ½ digit LED
Range	pH	0 to 14.00
	mV	0 to ±1999
Resolution	pH	0.01
	mV	1
Accuracy	pH	±0.01
	mV	±1
Temperature	–	Manual and Automatic (ATC)
Compensation	–	0 to 100°
Slope Correction	–	80 to 120%
Power	–	230v ± 10% ac, 50 Hz
Dimensions	–	275 x 175 x 76 mm



info@spectronicsindia.in



www.spectronicsindia.com



+91 9417350423

