



SPECTRONICS INDIA

AN ISO 9001:2015
CERTIFIED COMPANY

SI-221

DIGITAL TURBIDITY METER

We are offering our superior range of turbidity meter for our clients based globally. These are manufactured with the help of our designers and engineers using high technology and state-of-the-art facilities. Also, we are able to custom design the range in accordance with the specifications of the clients.



SALIENT FEATURES



HIGHLY STABLE & ACCURATE



3 1/2 DIGIT LED DISPLAY



30mm FLAT BOTTOM TEST TUBE



0 TO 1000 NTU

ACCESSORIES

- Test tubes
- Operation Manual
- Dust Cover

MODEL	SI-221 TABLE MODEL
Display	3 1/2 Digital LED
Detector	Photocell / Photodiode
Range	0 to 1000 NTU
Resolution	1 NTU
Accuracy	±3% FS± 1digit
Sample System	30 mm Clear Glass Test
Light Source	6V Tungsten Lamp or Bright Blue LED
Power	220V±10%, 50Hz AC
Dimensions	230mm x 180mm x 110mm



spectronicsindia@gmail.com



spectronicsindia.in



+91 9417350423

