



SPECTRONICS INDIA

**AN ISO 9001:2015
CERTIFIED COMPANY**

SI-223

DIGITAL TURBIDITY METER

We are offering our superior range of Turbidity Meter (with Nephelometer range) for our clients based globally. These are manufactured with the help of our designers and engineers using high technology and state-of-the-art facilities. Also, we are able to custom design the range in accordance with the specifications of the clients.



SALIENT FEATURES

- Highly Stable & Accurate
- 3½ Digit LED Display
- 30 mm flat bottom Test Tube
- Range 0 to 1000 NTU
(Selectable 0-200, 200-1000)

ACCESSORIES

- Test tubes
- Operation Manual
- Dust Cover

MODEL	SI-223
Type	Table Model
Display	3½ Digital LED
Detector	Photocell / Photodiode
Range	0 to 1000 NTU (With 0 to 200 and 200-1000 NTU)
Resolution	0.01 (0 to 200) & 1 NTU (200 to 1000)
Accuracy	±2% FS± 1digit
Sample System	30 mm Clear Glass Test
Light Source	6V Tungsten Lamp or Bright Blue LED
Power	220V±10%, 50Hz AC
Dimensions	230mm x 180mm x 110mm