

DIGITAL VISCOMETER

(TOUCH SCREEN)



SI- V205

Product Description

Digital Viscometer is a touch screen compact and easy to operate instruments with auto range showing time function for measurements. This is a unique system that combines convenience with accuracy and precision in analysis.

Features

Touch screen compact and easy to operate
Auto range showing time function for measurements
Sound alarm at under 20% torque, linear calibration

Applications

- Oil Paints
- Inks Latex Adhesives
- Juices and Dairy Products
- Pharmaceuticals
- Coating Solvents
- Polymer Solution

Technical Specification

Model	SI - V205
Measuring Range m.P.A.S.	10 - 2,000,000 mPAS
RPM	0.3, 0.6, 1.5, 3, 6, 12, 30, 60
No. of Standard Spindles	#1, #2, #3, #4 is standard configuration (#0 is for option)
Measurement Accuracy	±10 % of range
Repeatability	±0.5 %
Display	LCD Touch Screen
Temperature RTD Monitoring Sensor Probe	Yes
Output	RS 232 Interface
Power Supply	AC 220V/50 Hz
Dimensions L X H X B (mm)	400 X 350 X 360
Weight	9 Kg

OUR PRODUCT RANGE

Double Beam UV/Vis Spectrophotometer, Single Beam UV/Vis Spectrophotometer, Flame Photometer, Colorimeter, pH Meter, Conductivity Meter, DO Meter, Salinity Meter, Turbidity Meter, Colony Counter, Auto Karl-Fischer, Friability Apparatus, Milling Point Apparatus, Temperature Indicator, Water & Soil Analysis Kit, Portable pH Meter, Student Microscope, Binocular Microscope, Medical Microscope and Portable Instruments.



SPECTRONICS INDIA

Manufacturer of Scientific and Laboratory Instruments