



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

HEATING MANTLE

2000ML

MODEL SI 475



Suitable for laboratories, chemicals, pharmaceuticals, Dyes, Glass plants etc.



Particularly useful for Heating inflammable liquids.



TECHNICAL SPECIFICATION

- Heating Mantle 2000ML in 220 Volt
- Suitable for laboratories, chemicals, pharmaceuticals, Dyes, Glass plants etc.
- Particularly useful for Heating inflammable liquids.
- Electric heating net is hand knitted from glass Yarn and there is no joint to maintain the strength of the yarn.
- Stands high temperature up to 350 degree c.
- The body of the mantle is spun one piece from non-rusting Aluminum (up to 2.0 liters Capacity) duly stoving painted.
- All heating mantles are lagged with special grade mineral wool.



Electric heating net is hand knitted from glass Yarn and there is no joint to maintain the strength of the yarn.



The body of the mantle is spun one piece from non-rusting Aluminum (up to 2.0 liters Capacity) duly stoving painted.



All heating mantles are lagged with special grade mineral wool.

PRODUCT RANGE

Double Beam UV/VIS Spectrophotometer • Single Beam UV/VIS Spectrophotometer • Flame Photometer • Colorimeter
pH Meter • Conductivity Meter • DO Meter • Salinity Meter • Turbidity Meter • Colony Counter • Auto Karl- Fischer
Friability Apparatus • Melting Point Apparatus • Temperature Indicator • Water & Soil Analysis Kit • Portable pH Meter
Student Microscope • Binocular Microscope • Medical Microscope • Portable Instruments



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY



✉ spectronicsindia@gmail.com

🌐 spectronicsindia.in

☎ +91 9417350423