



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

HOT AIR OVEN

MODEL SI-445

With Double Wall construction, Inner made of Aluminum or Stainless steel and outer M.S. with powder coating with glass wool insulation, Temp. upto 250°C thermostatically controlled with built in L Shape thermometer. Supplied with adjustable shelves to on works on 220/230v AC.



SALIENT FEATURES

- Highly Stable & Accurate
- Electrical Operated
- Complete with three step Rotary Switch, Plug & cord.
- Temp. up to 250°C Shown in Top with Thermometer.
- Thermostat Controlled with Energy Regulator.
- Load & Main Indicator

SPECIFICATION

Temperature Range	Ambient to 250°C
Temperature Control	Thermostat Controlled with Energy Regulator
Temperature Indicator	Built in L Shape Thermometer
Inner Chamber	Aluminum or Stainless Steel
Outer Body	Mild Steel with Powder Coating
Insulation	Glass Wool
Shelves	2 Nos. (Adjustable)
Capacity	28 Ltrs.
Internal Size (HWD)	300 x 300 x 300 mm (12 x 12 x 12 inch)
Power Supply	220/230V AC, 50 Hz
Load	1 KV



Double Wall Construction



Temp. up to 250°C



Glass Wool Insulation



Thermostat Controlled



Energy Efficient

PRODUCT RANGE

- Double Beam UV/VIS Spectrophotometer • Single Beam UV/VIS Spectrophotometer • Flame Photometer • Colorimeter
- pH Meter • Conductivity Meter • DO Meter • Salinity Meter • Turbidity Meter • Colony Counter • Auto Karl- Fischer
- Friability Apparatus • Melting Point Apparatus • Temperature Indicator • Water & Soil Analysis Kit • Portable pH Meter
- Student Microscope • Binocular Microscope • Medical Microscope • Portable Instruments



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

