



**SPECTRONICS INDIA**

AN ISO 9001:2015 & CE CERTIFIED COMPANY

# DIGITAL HOT AIR OVEN

**MODEL SI-446**

With Double Wall construction, Inner made of Aluminum or Stainless steel and outer M.S. with powder coating with glass wool insulation, Temp. upto 250°C Display with Digital PID Controller. Supplied with adjustable shelves to on works on 220/230v AC.



## SALIENT FEATURES

- Highly Stable & Accurate
- Electrical Operated
- Complete with three step Rotary Switch, Plug & cord.
- Temp. up to 250°C Display with Digital PID Controller
- Thermostat Controlled with Energy Regulator.
- Load & Main Indicator

## SPECIFICATION

Temperature Range	Ambient to 250°C
Temperature Control	Digital PID Controller with Energy Regulator
Temperature Display	Digital
Inner Chamber	Aluminum or Stainless Steel
Outer Body	Mild Steel with Powder Coating
Insulation	Glass Wool
Shelves	2 Nos. (Adjustable)
Capacity	28 Ltrs.
Internal Size (HWD)	300 x 300 x 300 mm (12 x 12 x 12 inch)
Power Supply	220/230V AC, 50 Hz
Load	1 KV



DOUBLE WALL  
CONSTRUCTION



TEMP. UP TO  
250°C



DIGITAL PID  
CONTROLLER



GLASS WOOL  
INSULATION



ACCURATE &  
STABLE



ENERGY  
EFFICIENT

## PRODUCT RANGE

Double Beam UV/VIS Spectrophotometer • Single Beam UV/VIS Spectrophotometer • Flame Photometer • Colorimeter  
 • pH Meter • Conductivity Meter • DO Meter • Salinity Meter • Turbidity Meter • Colony Counter • Auto Karl- Fischer  
 • Friability Apparatus • Melting Point Apparatus • Temperature Indicator • Water & Soil Analysis Kit • Portable pH Meter  
 • Student Microscope • Binocular Microscope • Medical Microscope • Portable Instruments



**SPECTRONICS INDIA**

AN ISO 9001:2015 & CE CERTIFIED COMPANY



spectronicsindia@gmail.com



spectronicsindia.in



+91 9417350423