



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

Digital Melting point apparatus

MODEL **SI 225**



SALIENT FEATURES

- 4 Digit LED Display
- Direct Read-out facility
- Temp. 2°C above RT 350°C
- Readability 1°C
- Facility to hold last reading to complete or Power failed.
- Oil Free Instrument

TECHNICAL SPECIFICATION

1	Electronic Controller for adjusting the heating rate
2	Digital Display of temperature with PT-100 Sensor
3	Provision for automatic holding the melting point display
4	Temp range 2°C above RT 350°C
5	Readability 1°C
6	Display 4 Digit LED
7	Temp Sensor PT 100
8	Temperature setting gain & high / low facility
9	Stop led & hold led working melt position
10	Ready to start button use for sample testing time sample Fill and start
11	RUN LED Start sample testing running
12	P V To check the current/process value during "HOLD" condition
13	FAN Testing sample melt than fan running
14	Dimension : 325x265x170mm
15	Weight : 4.2Kg approx.
16	Power : 230V±10% AC 50 Hz

* We known as for Quality & Assured delivery to our valuable customer's & Users of our products. *

PRODUCT RANGE

Double Beam UV/VIS Spectrophotometer • Single Beam UV/VIS Spectrophotometer • Flame Photometer • Colorimeter
pH Meter • Conductivity Meter • DO Meter • Salinity Meter • Turbidity Meter • Colony Counter • Auto Karl- Fischer
Friability Apparatus • Melting Point Apparatus • Temperature Indicator • Water & Soil Analysis Kit • Portable pH Meter
Student Microscope • Binocular Microscope • Medical Microscope • Portable Instruments



SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY



spectronicsindia@gmail.com



spectronicsindia.in



+91 9417350423