



SPECTRONICS INDIA

AN ISO 9001:2015
CERTIFIED COMPANY

SI-167 DIGITAL COLONY COUNTER

Microprocessor Based Digital Colony Counter is a suitable instrument for counting Bacterial colonies. Our diverse range of Microprocessor Colony Counter offered by us is of premium quality and is fabricated using quality raw material. With the assistance of our designers and experts we are capable to offer customized solutions to our clients as per their specifications in stipulated time frame.



SALIENT FEATURES

-  Microprocessor Based counting control by Soft Touch key.
-  Audible Confirmation of each count
-  110mm Dia Magnifying Glass
-  Uniform glare free illumination
-  Digital Readout (4 Digits) from 0 – 999999
-  Auto Marker pen
-  Wolffhuegel glass grid with focusing facility

ACCESSORIES

- Marking Pen
- Magnifier Lens
- Dust Cover
- Manual

TECHNICAL SPECIFICATIONS

| PARAMETER | MODEL SI-167 |
|------------------|---|
| DIGITAL DISPLAY | 6 digit, 0 - 999999 Maximum Count |
| SOFT KEY CONTROL | Microprocessor based with facility of Counting control Increase/decrease delition |
| DISH SIZE | 110 mm |
| MAGNIFICATION | X 1.7 |
| DIMENSIONS | L 265 X B 230 X H 130 mm (Approx.) |
| WEIGHT | 2 Kg. (Approx.) |
| POWER | 230V ± 10% AC, 50Hz, 40W |
| ACCESSORIES | Marking Pen – (1 No.), Magnifier Lens – (1 Set) |