




SPECTRONICS INDIA

AN ISO 9001:2015 &
CE CERTIFIED COMPANY

SI-191

PORTABLE CONDUCTIVITY METER

We offer a wide range of Portable Conductivity Meter that is fabricated using quality raw material and is available at optimum prices. With the assistance of our engineers and designers we are able to carve a niche for ourselves in catering conductivity meter.



SALIENT FEATURES



Reliable and easy to operate



3 ½ Digit LCD Display



Cell Constant Adjustment Facility

ACCESSORIES

- CONDUCTIVITY CELL
- OPERATION MANUAL



CONDUCTIVITY CELL



TECHNICAL SPECIFICATIONS

PARAMETER	MODEL SI-191
Model	SI-191 Portable Model
Display	3½ Digit LCD
Range	0-1999µS 0-19.99 mS
Cell Constant Adjustment	Pre Adjustable
Conductivity Cell	Platinum DIP Type
RESOLUTION	0.1µS
Accuracy	±1% FS ±1 Digit
Power	1 x 9V Battery
Dimensions	165 x 85 x 28 mm
Weight	0.8 kg. (Approx.)
Accessories	Conductivity Cell, Operation Manual.



info@spectronicsindia.in



Spectronicsindia.com



+91 9417350423

