



SPECTRONICS INDIA

AN ISO 9001:2015 &
CE CERTIFIED COMPANY

SI-136

DIGITAL pH METER

The diverse range of Digital pH Meter catered by us is of top-notch quality and is manufactured using quality raw material and advanced technology. Our experts and engineers are capable enough to maintain pace with the latest trend of the marketplace. Further, we can custom design the range as per the specifications of the clients in the stipulated time frame.



Portable pH Meter SI-136



Combined pH Electrode

ACCESSORIES

- Combination pH Electrode,
- Buffer Tablets 4pH & 7pH



SALIENT FEATURES

- Highly Stable & Accurate
- LCD Display
- With 9V Battery

TECHNICAL SPECIFICATIONS

PARAMETER		SI-136 PORTABLE MODEL
Display		LCD Display
Range	pH	0 to 14.00
	mV	-
Resolution	pH	±0.01
	mV	-
Accuracy	pH	±0.01
	mV	-
Compensation		Manual 0 to 100° C
Power		1 × 9 V Battery
Dimensions		165 × 85 × 28 mm