




# SPECTRONICS INDIA

AN ISO 9001:2015 & CE CERTIFIED COMPANY

## PORTABLE TDS METER

### SI-186

We offer a wide range of Portable TDS Meter that is fabricated using quality raw material and is available at optimum prices. With the assistance of our engineers and designers we are able to carve a niche for ourselves in catering TDS meter.



#### SALIENT FEATURES



Reliable and easy to operate



LCD Display



Compact ABS Body

#### ACCESSORIES

- TDS CELL
- OPERATION MANUAL



Compact & Portable



Low Power Consumption



Accurate Results

#### SPECIFICATIONS

Model	SI-186 Portable Model
Display	3½ Digit LCD
Range	0-1999ppm (mg/L) 0-19.99 ppt (mg/L)
TDS Cell	Platinum DIP Type
RESOLUTION	0.1ppm
Accuracy	±1% FS ±1 Digit
Power	1 x 9V Battery
Dimensions	165 x 85 x 28 mm
Weight	0.8 kg. (Approx.)
Accessories	TDS Cell, Operation Manual.

info@spectronicsindia.in

www.spectronicsindia.com

+91 9417350423

